

TABLE OF CONTENTS

| | |
|-----------|--------------------------------|
| 03 | AMBULATORY INFUSION PUMP |
| 04 | ADMINISTRATION SETS |
| 05 | ACCESSORIES |
| 06 | ORDERING & CONTACT INFORMATION |

AMBULATORY INFUSION PUMP **NIMBUS™ II EMS**


REVOLUTIONARY ELECTRONIC INFUSION PUMP FOR **EMS.**

The Nimbus™ II EMS is a small, lightweight, smart, and durable infusion pump, making it a perfect candidate for emergency transport use.

This pump is highly customizable with editable presets, protocols and global safeguards, and names. Included are industry standard operating modes such as weight based programming and continuous volume over time.

IT1078

NIMBUS™ II EMS





| PRODUCT ID | PRODUCT NAME | QTY |
|------------|--|-----|
| IT1078 | NIMBUS™ II EMS PUMP INCLUDES: CLINCIAN IFU MANUAL | 1 |
| SOL222-2 | Preventative Maintenance | 1 |

NIMBUS™ II EMS ADMINISTRATION SETS

HS-005-B

QTY: 20/CASE





Tip protector, One Luer lock cap, One Female Luer Lock, One Cassette, One Anti-Free-Flow Clamp, One 0.2 Micron Air Eliminating Baby Filter, One Slide Clamp, One Backcheck Valve, One Male Luer Lock Adaptor, and Luer Cap.

Length: 67" (170cm)

HS-006

QTY: 20/CASE



Spike cap, One Spike, One Cassette, One Anti-Free-Flow Clamp, One 1.2 Micron Air Eliminating Filter, One Slide Clamp, One Backcheck Valve, One Male Luer Lock Adaptor, and Luer Cap.

Length: 90.5" (230 cm)

The Nimbus Administration Set is designed to administer fluids/medication from a container to a patient through a needle or catheter. The fluid path is sterile (EtO sterilized), non-pyrogenic, not made of DEHP, not made with natural latex. The set is single use only.

NIMBUS™ II EMS ACCESSORIES

CARRY POUCHES

QTY: EACH



SOL226

Fits Nimbus Pump with **500mL infusion bag**

OTHERS



IT1069

Nimbus Battery
EVE ER18505M 3.6V lithium battery

QTY: EACH or 24/SLEEVE



SOL234

Reinforced IV Pole Mount

QTY: EACH



SOL233

Single pump IV Pole Mount

QTY: EACH

INFORMATION

HOW TO ORDER & CONTACT US

OUR CUSTOMER SERVICE TEAM

Phone +1 (508) 650-2008
Fax +1 (508) 650-2006
Email orders@zynosolutions.com
Web www.zynosolutions.com

Operating days and hours: Monday - Friday, 8 AM - 5 PM

OUR TECH SUPPORT TEAM

Phone +1 (866) 395-1888
Email techsupport@zynosolutions.com
Web www.zynosolutions.com

Operating days and hours: 24/7

PREVENTATIVE MAINTENANCE

Customer is responsible for annual preventive maintenance procedures as required by the manufacturer. Zyno Solutions cannot assure safe use of the infusion pump if the customer fails to follow such procedures.

ADMINISTRATION SETS

Only manufacturer tubing is tested and calibrated for its products.



+1 (508) 650-2008
info@zynosolutions.com
177 Pine Street Natick MA 01760 USA
www.zynosolutions.com

ZYS-CT2E-003 Rev B