

**Z-800F**  
**Z-800W**  
**Z-800WF**

## IV Administration Set Catalog

Z-800F, Z-800W and Z-800WF

Zyno Administration Set Naming Convention .....3

B2/BX Primary Administration Sets: Z-800F/W/WF.....4

B2/BX Primary Administration Filter Sets: Z-800F/W/WF.....5

Secondary Administration Sets.....6

Extension Administration Sets with Filter.....7

Administration Set Quick Reference for: Z-800F, Z-800W, Z-800WF.....8

Ordering Information.....9


## Z-800 Series Large Volume Infusion Pumps Administration Sets

These administration sets support the following Large Volume Infusion Pumps:

- Zyno Z-800F Set Based Anti-free Flow Pump
- ZynoFlo Z-800W Set Based Anti-free Flow Pump
- ZynoFlo Z-800WF Set Based Anti-free Flow Pump

## Administration Set Naming Convention

The naming convention allows for the recognition of the administration set via the part number. The naming convention will identify the OEM manufacturer, the pump type, with or without needleless injection site(s), location of injection site and filter options.











**Warning:** Only use Zyno Medical administration sets provided by Zyno Solutions with the Z-800 Series Large Volume Infusion Pumps. There are risks associated with using any administration sets other than Zyno Medical administration sets with its device. Zyno Medical's warranty for its device will be null and void and Zyno Medical will assume no responsibility for any incidents that may occur if the device is not utilized strictly in accordance with its product labeling.

# Primary B2/BX Administration Sets

## B2/BX Primary Administration Sets: Z-800F/W/WF

### B2/BX Primary Administration Sets

B2/BX Primary Administration Sets are single-use sets for use with Zyno Medical Z-800F, Z-800W and Z-800WF Infusion pumps. All B2/BX Primary administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.



Product #	Description	Length	# of Y-Sites	Qty
B2-70070 or BX-70070	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Pinch clamp, (5)-Roller clamp, (6)-Male Luer lock adaptor, (7)-Luer lock tip protector.</p> <p>Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p> 	105" (266cm)	0	50
B2-70071-D or BX-70071-D	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Pinch clamp, (5)-Roller clamp, (6)-Slide clamp, (7)-Needleless injection site, (8)-Male Luer lock adaptor, (9)-Luer lock tip protector.</p> <p>Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p> 	105" (266cm)	1 (distal)	50
B2-70071-P or BX-70071-P	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5)-Needleless injection site, (6)-Pinch clamp, (7)-Roller clamp, (8)-Male Luer lock adaptor, (9)-Luer lock tip protector.</p> <p>Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p> 	105" (266cm)	1 (proximal)	50
B2-70072 or BX-70072	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5 &amp; 9)-Needleless injection sites, (6)-Pinch clamp, (7)-Roller clamp, (8)-Slide clamp, (10)-Male Luer lock adaptor, (11)-Luer lock tip protector.</p> <p>Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p> 	105" (266cm)	2	50

# Primary B2/BX Administration Sets

## B2/BX Primary Administration Filter Sets: Z-800F/W/WF


### B2/BX Primary Sets with 0.2 Micron Filter

B2/BX Primary Administration Sets are single-use sets for use with Zyno Medical Z-800F, Z-800W and Z-800WF Infusion pumps. All B2/BX Primary administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.

Product #	Description	Length	# of Y-Sites	Qty
B2-70071-DF or BX-70071-DF	 (1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Pinch clamp, (5)-Roller clamp, (6)-0.2 micron filter, (7)-Slide clamp, (8)-Needleless injection site, (9)-Male Luer lock adaptor, (10)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL	105" (266cm)	0	50
B2-70072-F or BX-70072-F	 (1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5 & 10)-Needleless injection sites, (6)-Pinch clamp, (7)-Roller clamp, (8)-0.2 micron filter, (9)-Slide clamp, (11)-Male Luer lock adaptor, (12)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL	105" (266cm)	1 (distal)	50

### B2/BX Primary Administration Sets with 1.2 Micron Filter

B2/BX Primary Administration Sets are for use with Zyno Medical Z-800F, Z-800W and Z-800WF Infusion pumps. All B2/BX Primary administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.



Product #	Description	Length	# of Y-Sites	Qty
B2-70071-DF120 or BX-70071-DF120	 (1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5 & 9)-Needleless injection sites, (6)-Pinch clamp, (7)-Roller clamp, (8)-Slide clamp, (10)-Male Luer lock adaptor, (11)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL	105" (266cm)	2	50

# Secondary Administration Sets

## Secondary Administration Sets

### Secondary Administration Sets

Secondary Administration Sets are single-use sets for use with Zyno Medical Z-800 Infusion pump. All Secondary administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.


Product #	Description	Length	# of Y-Sites	Qty
AA-80075	 (1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-Male Luer lock adaptor, (6)-Luer lock tip protector, (7)-Bag hanger (11 inch). Priming Volume: 9 mL (approx.)/ Drop Size: 20 drops/mL	40" (101cm)	0	50
AX-80075	 (1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-Male Luer lock adaptor, (6)-Luer lock tip protector, (7)-Bag hanger (16.5 inch). Priming Volume: 9 mL (approx.)/ Drop Size: 20 drops/mL	40" (101cm)	0	50

# Primary B2/BX Administration Sets

## Extension Administration Sets with Filter


### Extension Administration Sets with 0.2 Micron Filter

Extension Administration Sets are single-use sets for use with Zyno Medical Z-800F, Z-800W and Z-800WF Infusion pumps. All extension administration set fluid paths are sterile, non-Pyrogenic are not made with natural rubber latex, and are not made with DEHP.

Product #	Description	Length	# of Y-Sites	Qty
FE01		10" (25cm)	0	50
(1)-Tip protector, (2)-0.2 micron filter with female Luer adapter, (3)-Slide clamp, (4)-Male Luer lock adaptor, (5)-Luer lock tip protector.				

### Extension Administration Sets with 1.2 Micron Filter

Extension Administration Sets are single-use sets for use with Zyno Medical Z-800F, Z-800W and Z-800WF Infusion pumps. All extension administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.

Product #	Description	Length	# of Y-Sites	Qty
FE04		14" (35cm)	0	50
(1)-Tip protector, (2)-Female Luer adapter, (3)-1.2 micron filter, (4)-Slide clamp, (5)-Male Luer lock adaptor, (6)-Luer lock tip protector.				

# B2/BX Administration Sets Quick Reference

Z-800 Series Large Volume Infusion Pumps Quick Reference for: Z-800F, Z-800W, Z-800WF												
Product Number	Primary Admin Sets	Secondary Admin Sets	Filter Extension	Back-check valve	Filter (μ)	Y-Site(s)		Roller Clamp	Slide Clamp	Drip chamber	Length	Units/ Case
						Proximal	Distal					
B2-70070 BX-70070	X							X		X	105"	50
B2-70071-D BX-70071-D	X						X	X	X	X	105"	50
B2-70071-P BX-70071-P	X			X		X		X		X	105"	50
B2-70071-DF BX-70071-DF	X				0.2		X	X	X	X	105"	50
B2-70072 BX-70072	X			X		X	X	X	X	X	105"	50
B2-70072-F BX-70072-F	X			X	0.2	X	X	X	X	X	105"	50
AA-80075 AX-80075		X						X		X	40"	50
FE01			X		0.2				X		10"	50
FE04			X		1.2				X		14"	50
B2-70071-DF120 BX-70071-DF120	X				1.2		X	X	X	X	105"	50
Administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.												



# Ordering Information

## Zyno Solutions

Zyno Solutions is a trusted provider of intelligent infusion systems for the alternate care market. We provide industry leading products where there is a need for durable and advanced technology infusion systems.

The technologies that we offer provide the highest quality of care and patient safety along with our superior service and support. Our goal is to provide infusion system designed specifically for the infusion needs of our customers.

## Preventative Maintenance

Customer is responsible for annual preventive maintenance procedures as required by the manufacturer. Zyno Solutions cannot assure safe use of the infusion pump if the customer fails to follow such procedures.

## Administration Sets

Only manufacturer tubing is tested and calibrated for its products.

## Customer Service Contact Information

Phone: 508-650-2008  
Fax: 508-650-2006  
[orders@zynosolutions.com](mailto:orders@zynosolutions.com)  
[www.zynosolutions.com](http://www.zynosolutions.com)

## Technical Support

866-395-1988  
[techsupport@zynosolutions.com](mailto:techsupport@zynosolutions.com)



Zyno Solutions, LLC  
177 Pine Street  
Natick, MA 01760  
Phone: 866-395-1988  
Email: [info@zynosolutions.com](mailto:info@zynosolutions.com)