

Z-800F/W/WF

# ADMINIS- TRATION SETS

---

## **2020 PRODUCT CATALOG**

Includes B2/BX primary, secondary and extension administration sets for Z-800F, Z-800W, and Z-800WF.

**P** 1 508 650 2008

**E** [orders@zynosolutions.com](mailto:orders@zynosolutions.com)  
**www.zynosolutions.com**



# TABLE OF CONTENT

---

About Zyno Solutions	03
Zyno Administration Set Naming Convention	04
Quick Reference Guide	05
B2/BX Primary Administration Sets	06
Secondary Administration Sets	08
Extension Administration Sets	08
Ordering Information	09



## WE ARE **ZYNO SOLUTIONS.**

Zyno Solutions is a trusted provider of intelligent infusion systems for the alternate care market. We provide industry leading products where there is a need for durable and advanced technology infusion systems.

The technologies that we offer provide the highest quality of care and patient safety along with our superior service and support. Our goal is to provide infusion systems designed specifically for the infusion needs of our customers.

Visit [www.zynosolutions.com](http://www.zynosolutions.com) for more information.

# ABOUT Z-800F/W/WF ADMINISTRATION SETS


This catalog is for Zyno Medical Large Volume Infusion Pumps:

- Z-800F set-based, anti-free flow pump
- ZynoFlo Z-800W pump
- ZynoFlo Z-800WF set-based, anti-free flow pump



## HOW ADMINISTRATION SETS ARE NAMED

The naming convention allows for the recognition of the administration set via the part number. The naming convention will identify the OEM manufacturer, the pump type, with or without needleless injection site(s), location of injection site and filter options.



**WARNING:** Only use Zyno Medical provided administration sets with the Z-800 Series Large Volume Infusion Pumps. There are risks associated with using any administration sets other than Zyno Medical administration sets with this device. Zyno's warranty for its device will be null and void and Zyno will assume no responsibility for any incidents that may occur if the device is not utilized strictly in accordance with its product labeling.

# QUICK REFERENCE FOR ADMINISTRATION SETS







Product Number	Primary Set	Secondary Set	Filter Extension	Back-check valve	Filter	Y-Site(s)		Roller Clamp	Slide Clamp	Drip chamber	Length	Units/Case
						proximal	distal					
B2-70070 BX-70070	X							X		X	105"	50
B2-70071-D BX-70071-D	X						X	X	X	X	105"	50
B2-70071-P BX-70071-P	X			X		X		X		X	105"	50
B2-70071-DF BX-70071-DF	X				0.2		X	X	X	X	105"	50
B2-70072 BX-70072	X			X		X	X	X	X	X	105"	50
B2-70072-F BX-70072-F	X			X	0.2	X	X	X	X	X	105"	50
AA-80075 AX-80075		X						X		X	40"	50
FE01			X		0.2				X		10"	50
FE04			X		1.2				X		14"	50
B2-70071-DF120 BX-70071-DF120	X				1.2		X	X	X	X	105"	50

- Administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP
- Single-use sets, universal spike, spike protector, drip chamber, male luer-lock and tip protector
- Needless injection sites
- Drip chamber = 20 drops/ mL

# ADMINISTRATION SETS


All of the administration listed below are single-use sets for use with Zyno Medical Z-800F, Z-800W and Z-800WF Infusion pumps. The administration sets (B2/BX primary administration sets, secondary administration sets, and extension administration sets) fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.

## B2/BX PRIMARY ADMINISTRATION SETS

PRODUCT ID	DESCRIPTION	LENGTH	# OF Y-SITES	QTY
B2-70070	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Pinch clamp, (5)-Roller clamp, (6)-Male Luer lock adaptor, (7)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	0	50
or				
BX-70070				
B2-70071-D	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Pinch clamp, (5)-Roller clamp, (6)-Slide clamp, (7)-Needleless injection site, (8)-Male Luer lock adaptor, (9)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	1 (distal)	50
or				
BX-70071-D				
B2-70071-P	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5)-Needleless injection site, (6)-Pinch clamp, (7)-Roller clamp, (8)-Male Luer lock adaptor, (9)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	1 (proximal)	50
or				
BX-70071-P				



B2-70072



or

(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5 & 9)-Needleless injection sites, (6)-Pinch clamp, (7)-Roller clamp, (8)-Slide clamp, (10)-Male Luer lock adaptor, (11)-Luer lock tip protector.


Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL

105"  
(266cm)

2

50

BX-70072




## B2/BX PRIMARY ADMINISTRATION SETS WITH 0.2 MICRON FILTER

PRODUCT ID DESCRIPTION

LENGTH # OF Y-SITES QTY

B2-70071-DF



or

(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Pinch clamp, (5)-Roller clamp, (6)-0.2 micron filter, (7)-Slide clamp, (8)-Needleless injection site, (9)-Male Luer lock adaptor, (10)-Luer lock tip protector.


Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL

105"  
(266cm)


1 (distal)

50

BX-70071-DF



B2-70072-F



or

(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5 & 10)-Needleless injection sites, (6)-Pinch clamp, (7)-Roller clamp, (8)-0.2 micron filter, (9)-Slide clamp, (11)-Male Luer lock adaptor, (12)-Luer lock tip protector.


Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL

105"  
(266cm)



2

50



BX-70072-F




## B2/BX PRIMARY ADMINISTRATION SETS WITH 1.2 MICRON FILTER

PRODUCT ID	DESCRIPTION	LENGTH	# OF Y-SITES	QTY
B2-70071-DF120	 <p>(1)- Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Pinch clamp, (5)-Roller clamp, (6)-1.2 micron filter, (7)-Slide clamp, (8)-Needleless injection site, (9)-Male Luer lock adaptor, (10)-Luer lock tip protector.</p>	105" (266cm)	1 (distal)	50
or	Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL			
BX-7071-DF120				

## SECONDARY ADMINISTRATION SETS


PRODUCT ID	DESCRIPTION	LENGTH	# OF Y-SITES	QTY
AA-80075	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-Male Luer lock adaptor, (6)-Luer lock tip protector, (7)-Bag hanger (11 inch).</p>	40" (101cm)	0	50
	Priming Volume: 9 mL (approx.)/ Drop Size: 20 drops/mL			
AX-80075	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-Male Luer lock adaptor, (6)-Luer lock tip protector, (7)-Bag hanger (16.5 inch).</p>	40" (101cm)	0	50
	Priming Volume: 9 mL (approx.)/ Drop Size: 20 drops/mL			

## EXTENSION ADMINISTRATION SETS WITH 0.2 MICRON FILTER

PRODUCT ID	DESCRIPTION	LENGTH	# OF Y-SITES	QTY
FE01	 <p>(1)-Tip protector, (2)-0.2 micron filter with female Luer adapter, (3)-Slide clamp, (4)-Male Luer lock adaptor, (5)-Luer lock tip protector.</p>	10" (25cm)	0	50



# EXTENSION ADMINISTRATION SETS WITH 1.2 MICRON FILTER

PRODUCT ID	DESCRIPTION	LENGTH	# OF Y-SITES	QTY
FE04	<div> <div> <div>1</div> <div>2</div> <div>3"</div> <div>3</div> <div>6"</div> <div>4</div> <div>5</div> <div>6</div> </div>  <p>(1)-Tip protector, (2)-Female Luer adapter, (3)-1.2 micron filter, (4)-Slide clamp, (5)-Male Luer lock adaptor, (6)-Luer lock tip protector.</p> </div>	14" (35cm)	0	50



Ready to Order? Need Help?  
Let’s get in touch.

### OUR CUSTOMER SERVICE TEAM

**P** 1 (508) 650 2008  
**F** 1 (508) 650-2006  
**E** [orders@zynosolutions.com](mailto:orders@zynosolutions.com)  
[www.zynosolutions.com](http://www.zynosolutions.com)

Hours: 8 AM - 5 PM, Monday to Friday

For technical support, please contact 866-395-1988 or [techsupport@zynosolutions.com](mailto:techsupport@zynosolutions.com).

