



# TABLE OF CONTENTS

<b>03</b>	AMBULATORY INFUSION PUMP
<b>04</b>	SPIKE SYSTEM
<b>05</b>	LUER LOCK SYSTEM
<b>08</b>	ACCESSORIES
<b>11</b>	ORDERING & CONTACT INFORMATION

# AMBULATORY INFUSION PUMP NIMBUS™ II FLEX/PLUS

## THE FUTURE OF AMBULATORY INFUSION THERAPY.


The Nimbus™ II Ambulatory Infusion System combines the accuracy of an electronic infusion pump with simple and intuitive programming.

This unique pump is engineered to be forward-thinking and capable of delivering advantages that meet the evolving needs of cancer care.




**IT1062**

NIMBUS™ II FLEX



**IT1078**

NIMBUS™ II PLUS




PRODUCT ID	PRODUCT NAME	QTY
IT1062	NIMBUS™ II FLEX PUMP	1
IT1078	NIMBUS™ II PLUS PUMP	1
<b>INCLUDED WITH EITHER PUMP:</b>		
SOL211	CLINICAL QUICK REFERENCE GUIDE	1
SOL212	PATIENT QUICK REFERENCE GUIDE	1
IT1069	NIMBUS BATTERY (EVE ER18505M 3.6V lithium battery)	1

# ADMINISTRATION SETS & IV BAGS

## SPIKE SYSTEM

**HS-004**

**QTY: 20/CASE**



(1) Luer Cap, (2) One Female Luer Lock Adaptor, (3) One Cassette, (4) One 1.2 Micron Filter, (5) One Slide Clamp, (6) One Backcheck Valve, (7) One Male Luer Lock Adaptor, (8) Luer Cap  
Priming Volume: 3mL (approx.)  
Length: 67" (170cm)

Our filters are made with low protein binding PES material.

**SECURE® IV BAGS**

**QTY: 50/CASE**



**SOL202**

50mL fill

**SOL201**

100mL fill

**SOL203**

250mL fill

**SOL206**

500mL fill

**SOL207**


1000mL fill

# ADMINISTRATION SETS & IV BAGS

## LUER LOCK SYSTEM

### HS-008



**QTY: 20/CASE**



(1) Luer Cap, (2) One Proximal Male Luer Lock Adaptor, (3) One Cassette, (4) One 0.2 Micron Filter, (5) One Slide Clamp, (6) One Backcheck Valve, (7) One Distal Male Luer Lock Adaptor, (8) Luer Cap  
Priming Volume: 2.5mL (approx.)  
Length: 53" (135cm)

### HS-008-B

**QTY: 20/CASE**



(1) Luer Cap, (2) One Proximal Male Luer Lock Adaptor, (3) One Cassette, (4) One 0.2 Micron Baby Filter, (5) One Slide Clamp, (6) One Backcheck Valve, (7) One Distal Male Luer Lock Adaptor, (8) Luer Cap  
Priming Volume: 2.5mL (approx.)  
Length: 53" (135cm)

Our filters are made with low protein binding **PES** material.

The Nimbus Administration Set is designed to administer fluids/medication from a container to a patient through a needle or catheter. The fluid path is sterile (EtO sterilized), non-pyrogenic, not made of DEHP, not made with natural latex. The set is single use only.

# ADMINISTRATION SETS & IV BAGS

## LUER LOCK SYSTEM

### EVA BAGS

QTY: 20/CASE



**SOL257**

150mL fill



**SOL258**

250mL fill



**SOL273**

500mL fill

Only SOL273 comes 50 bags/case.



**SOL274**

105mL fill (for mini hardshell)

# ADMINISTRATION SETS & IV BAGS

## LUER LOCK SYSTEM

**EVA BAGS 2 WAY NEEDLE-FREE VALVE** QTY: 20/CASE



**SOL260**

150mL fill



**SOL261**

250mL fill



**SOL262**

500mL fill



# CARRY POUCHES, HARDSHELL & MORE **ACCESSORIES**

## CARRY POUCHES

QTY: EACH



**SOL224**

Fits Nimbus Pump with  
**50mL Spike System** OR  
**100mL LL System**



**SOL247**

Fits Nimbus Pump with **100mL Spike System** OR **250mL LL System**  
**\*\*\*with belt loop\*\*\***



while  
supplies  
last

**SOL225**

Fits Nimbus Pump with **100mL Spike System** OR **250mL LL System**



**SOL226**

Fits Nimbus Pump with  
**250mL Spike System** OR  
**500mL LL System**



**SOL248**

Fits Nimbus Pump with **500mL Spike System** OR **1L LL System**



# CARRY POUCHES, HARDSHELLS & MORE **ACCESSORIES**

## **HARDSHELLS**

**QTY: 20/CASE**



**SOL228**

For 100mL size  
\*\*\*Max Bag fill is 110mL\*\*\*



**SOL244**

For 250mL size  
\*\*\*Max Bag fill is 250mL\*\*\*



**SOL254**

For mini EVA bag  
\*\*\*Max Bag fill is 105 mL \*\*\*

## **HARDSHELL POUCHES**

**QTY: EACH**



**SOL251**

For 100mL hardshell with  
window



**SOL252**

For 250mL hardshell with window



**SOL265**

For mini hardshell with window

# CARRY POUCHES, HARDSHELLS & MORE **ACCESSORIES**

## OTHERS



**IT1069**

**Nimbus Battery**

EVE ER18505M 3.6V lithium  
battery

**QTY: EACH or 24/SLEEVE**



**SOL234**

**Chemo Spill Kit** includes  
absorbant Bio-wipe Bag, Pr Chemo  
Gloves, Chemo Bag, Instructions

**QTY: 5/PACK**

# INFORMATION

# HOW TO ORDER & CONTACT US

## OUR CUSTOMER SERVICE TEAM

**Phone** +1 (508) 650-2008  
**Fax** +1 (508) 650-2006  
**Email** [flexorders@zynosolutions.com](mailto:flexorders@zynosolutions.com)  
**Web** [www.zynosolutions.com](http://www.zynosolutions.com)

**Operating days and hours:** Monday - Friday, 8 AM - 5 PM

## OUR TECH SUPPORT TEAM

**Phone** +1 (866) 395-1888  
**Email** [techsupport@zynosolutions.com](mailto:techsupport@zynosolutions.com)  
**Web** [www.zynosolutions.com](http://www.zynosolutions.com)

**Operating days and hours:** 24/7

## PREVENTATIVE MAINTENANCE

Customer is responsible for annual preventive maintenance procedures as required by the manufacturer. Zyno Solutions cannot assure safe use of the infusion pump if the customer fails to follow such procedures.

## ADMINISTRATION SETS

Only manufacturer tubing is tested and calibrated for its products.



+1 (508) 650-2008  
info@zynosolutions.com  
177 Pine Street Natick MA 01760 USA  
www.zynosolutions.com

ZYS-CT2F-002 Rev H