



**SIMPLE.
SAFE.
FLEXIBLE.**

Nimbus™ II Flex



WHY NIMBUS™?

The Nimbus Ambulatory Infusion System combines the accuracy of an electronic infusion pump with a simple and easy to understand user interface. Nimbus is uniquely engineered to be forward-thinking and capable of delivering advantages that meet the evolving needs of cancer care. Nimbus provides a system that uniquely meets current needs while evolving toward the future to help improve the patient outcomes.

NIMBUS FEATURES

- ✓ Electronically Programmed rate and volume delivery
- ✓ Durable for repeated ambulatory use
- ✓ Preloaded chemotherapy protocols
 - Practice standard protocols
 - General hours preloaded protocols
- ✓ Automatically presents entered Rx for review
- ✓ Ideal for low volume infusions
- ✓ Accuracy of +/-5%
- ✓ Meets CMS and DME standard



PORTABLE

The extremely compact design (6.14 ounces) provides a discrete profile for ambulatory patients.



SAFE

Integrated "Watchdog" software and hardware guards may help to prevent programming errors.



AFFORDABLE

Nimbus is exceedingly cost effective to maximize your ROI while also increasing efficiency.

Optimize the Cancer Care Workflow

Schedule a product demo



Get in touch:
info@zynosolutions.com

Follow Us:

www.zynosolutions.com



TECHNICAL SPECIFICATIONS

Pump Mechanism	Linear Peristaltic	
Operating Modes	Continuous Mode	
Maximum Infusion Pressure	30 psi	
Maximum Time to Occlusion Alarm	20 minutes	
Rate	0.1mL/h to 135 mL/h in 0.1 mL/h increments	
VTBI	1-1500 mL by 0.1 mL increment; single infusion VTBI maximum value limited by remaining pump lifetime VTBI	
Flow Rate Accuracy	+/- 5%	
Dimensions	4.25 in. x 2.28 in. x 1.57 in. (108 mm x 58 mm x 40 mm)	
Weight	6.14 ounces (174 grams)	
Power Supply	AC/DC power via Micro-USB port, internally powered by battery (battery not rechargeable)	
Battery life	400 h or 1500 mL of volume infused, whichever is reached first	
Alarms	<ul style="list-style-type: none"> • Cassettes Loading Error • System Error • Firmware Error • Upstream Occlusion • Downstream Occlusion • Battery Depleted 	<ul style="list-style-type: none"> • Invalid Infusion Parameter • Pump Unattended • Infusion Complete • Power off and replace battery • Due for Service
Warnings	<ul style="list-style-type: none"> • Battery Low 	
Administration Set	<ul style="list-style-type: none"> • Bag Spike • Free Flow Protection • Slide Clamp 	<ul style="list-style-type: none"> • 1.2 Micron Filter (air eliminating) • Luer Lock Adapter
Accessories	Carrying pouch for pump, drug bag, and hardshell	
Standards	Electrical Safety: IEC 60601-1 Pump Performance: IEC 60601-2-24 Electromagnetic Compatibility: IEC 60601-1-2 Usability and Human Factor: IEC 60601-1-6 ed3 / IEC 62366 Alarm System: IEC 60601-1-8 Home Healthcare Environment: IEC 60601-1-11	