

Z-800

IV Administration Set Catalog

Z-800

Ordering Information.....3

Zyno Administration Set Naming Convention.....4

A2/AX Primary Administration Sets: Z-800.....5

A2/AX Primary Filter Administration Sets: Z-800.....6-7

Secondary Administration Sets.....8

Administration Set Quick Reference for: Z-800.....9


Z-800 Series Large Volume Infusion Pumps Administration Sets

These administration sets support the following Large Volume Infusion Pumps:

➤ Zyno Z-800 Basic Infusion Pump

Administration Set Naming Convention

The naming convention allows for the recognition of the administration set via the part number. The naming convention will identify the OEM manufacturer, the pump type, with or without needless injection site(s), location of injection site and filter options.







Warning: Only use Zyno Medical administration sets provided by Zyno Solutions with the Z-800 Series Large Volume Infusion Pumps. There are risks associated with using any administration sets other than Zyno Medical administration sets with its device. Zyno Medical's warranty for its device will be null and void and Zyno Medical will assume no responsibility for any incidents that may occur if the device is not utilized strictly in accordance with its product labeling.

Primary A2/AX Administration Sets

A2/AX Primary Administration Sets: Z-800

A2/AX Primary Administration Sets

A2/AX Primary Administration Sets are single-use sets for use with Zyno Medical Z-800 Infusion pump. All A2/AX Primary administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.




Product #	Description	Length	# of Y-Sites	Qty
A2-80070 or AX-80070	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5) -Male Luer lock adaptor, (6)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	0	50
A2-80071-D or AX-80071-D	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-Slide clamp, (6)-Needleless injection site, (7)-Male Luer lock adaptor, (8)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	1 (distal)	50
A2-80071-P or AX-80071-P	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5)-Needleless injection site, (6)-Roller clamp, (7)-Male Luer lock adaptor, (8)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	1 (proximal)	50
A2-80072 or AX-80072	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5 & 8)-Needleless injection sites, (6)-Roller clamp, (7)-Slide clamp, (9)-Male Luer lock adaptor, (10)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	2	50

Primary A2/AX Administration Sets

A2/AX Primary Filter Administration Sets: Z-800



A2/AX Primary Administration Sets with 0.2 Micron Filter

A2/AX Primary Administration Sets are single-use sets for use with Zyno Medical Z-800 Infusion pumps. All A2/AX Primary administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.

Product #	Description	Length	# of Y-Sites	Qty
A2-80071-DF or AX-80071-DF	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-0.2 micron filter, (6)-Slide clamp, (7)-Needleless injection site, (8)-Male Luer lock adaptor, (9)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	1 (distal)	50
AX-80071-PF	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5)-Needleless injection site, (6)-Roller clamp, (7)-0.2 micron filter, (8)-Male Luer lock adaptor, (9)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	1 (proximal)	50
A2-80072-F or AX-80072-F	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5 & 9)-Needleless injection sites, (6)-Roller clamp, (7)-0.2 micron filter, (8)-Slide clamp, (10)-Male Luer lock adaptor, (11)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	2	50

A2/AX Primary Administration Sets with 1.2 Micron Filter



A2/AX Primary Administration Sets are single use sets for use with Zyno Medical Z-800 Infusion pumps. All A2/AX Primary administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.

Product #	Description	Length	# of Y-Sites	Qty
A2-80071-DF120	 (1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-1.2 micron filter, (6)-Slide clamp, (7)-Needleless injection site, (8)-Male Luer lock adaptor, (9)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL	105" (266cm)	1 (distal)	50
or				
AX-80071-DF120				

Secondary Administration Sets

Secondary Administration Sets

Secondary Administration Sets are single-use sets for use with Zyno Medical Z-800 Infusion pump. All Secondary administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.

Product #	Description	Length	# of Y-Sites	Qty
AA-80075		40" (101cm)	0	50
	(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-Male Luer lock adaptor, (6)-Luer lock tip protector, (7)-Bag hanger (11 inch).			
	Priming Volume: 9 mL (approx.)/ Drop Size: 20 drops/mL			
AX-80075		40" (101cm)	0	50
	(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-Male Luer lock adaptor, (6)-Luer lock tip protector, (7)-Bag hanger (16.5 inch).			
	Priming Volume: 9 mL (approx.)/ Drop Size: 20 drops/mL			

A2/AX Administration Sets Quick Reference

Z-800 Series Large Volume Infusion Pumps

Quick Reference for: Z-800

Product Number	Primary Admin Sets	Secondary Admin Sets	Filter Extension	Back-check valve	Filter (μ)	Y-Site(s)		Roller Clamp	Silde Clamp	Drip chamber	Length	Units/ Case
						Proximal	Distal					
A2-80070 AX-80070	X							X		X	105"	50
A2-80071-D AX-80071-D	X						X	X	X	X	105"	50
A2-80071-P AX-80071-P	X			X		X		X		X	105"	50
A2-80071-DF AX-80071-DF	X				0.2		X	X	X	X	105"	50
AX-80071-PF	X			X	0.2	X		X		X	105"	50
A2-80072 AX-80072	X			X		X	X	X	X	X	105"	50
A2-80072-F AX-80072-F	X			X	0.2	X	X	X	X	X	105"	50
AA-80075 AX-80075		X						X		X	40"	50
FE01			X		0.2				X		10"	50
FE04			X		1.2				X		14"	50
A2-80071-DF120 AX-80071-DF120	X				1.2		X	X	X	X	105"	50

Administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.

Ordering Information

Zyno Solutions

Zyno Solutions is a trusted provider of intelligent infusion systems for the alternate care market. We provide industry leading products where there is a need for durable and advanced technology infusion systems.

The technologies that we offer provide the highest quality of care and patient safety along with our superior service and support. Our goal is to provide infusion system designed specifically for the infusion needs of our customers.

Preventative Maintenance

Customer is responsible for annual preventive maintenance procedures as required by the manufacturer. Zyno Solutions cannot assure safe use of the infusion pump if the customer fails to follow such procedures.

Administration Sets

Only manufacturer tubing is tested and calibrated for its products.

Customer Service Contact Information

Phone: 508-650-2008
Fax: 508-650-2006
orders@zynosolutions.com
www.zynosolutions.com

Technical Support

866-395-1988
techsupport@zynosolutions.com



Zyno Solutions, LLC
177 Pine Street
Natick, MA 01760
Phone: 866-395-1988
Email: info@zynosolutions.com