

Z-800

ADMINIS- TRATION SETS

2020 PRODUCT CATALOG

Includes A2/AX primary,
secondary and extension
administration sets for Z-800.

P 1 508 650 2008

E orders@zynosolutions.com

www.zynosolutions.com

TABLE OF CONTENTS

About Zyno Solutions	03
Zyno Administration Set Naming Convention	04
Quick Reference Guide	05
A2/AX Primary Administration Sets	06
Secondary Administration Sets	08
Ordering Information	09



WE ARE **ZYNO SOLUTIONS.**

Zyno Solutions is a trusted provider of intelligent infusion systems for the alternate care market. We provide industry leading products where there is a need for durable and advanced technology infusion systems.

The technologies that we offer provide the highest quality of care and patient safety along with our superior service and support. Our goal is to provide infusion systems designed specifically for the infusion needs of our customers.

Visit www.zynosolutions.com for more information.

ABOUT Z-800 ADMINISTRATION SETS


This catalog is for Zyno Medical Large Volume Infusion Pumps:

- Zyno Z-800 Basic Infusion Pump



HOW ADMINISTRATION SETS ARE NAMED

The naming convention allows for the recognition of the administration set via the part number. The naming convention will identify the OEM manufacturer, the pump type, with or without needless injection site(s), location of injection site and filter options.



WARNING: Only use Zyno Medical provided administration sets with the Z-800 Series Large Volume Infusion Pumps. There are risks associated with using any administration sets other than Zyno Medical administration sets with this device. Zyno's warranty for its device will be null and void and Zyno will assume no responsibility for any incidents that may occur if the device is not utilized strictly in accordance with its product labeling.

QUICK REFERENCE FOR ADMINISTRATION SETS







Product Number	Primary Set	Secondary Set	Filter Extension	Back-check valve	Filter	Y-Site(s)		Roller Clamp	Slide Clamp	Drip chamber	Length	Units/Case
						proximal	distal					
A2-80070 AX-80070	X							X		X	105"	50
A2-80071-D AX-80071-D	X						X	X	X	X	105"	50
A2-80071-P AX-80071-P	X			X		X		X		X	105"	50
A2-80071-DF AX-80071-DF	X				0.2		X	X	X	X	105"	50
AX-80071-PF	X			X	0.2	X		X		X	105"	50
A2-80072 AX-80072	X			X		X	X	X	X	X	105"	50
A2-80072-F AX-80072-F	X			X	0.2	X	X	X	X	X	105"	50
AA-80075 AX-80075		X						X		X	40"	50
FE01			X		0.2				X		10"	50
FE04			X		1.2				X		14"	50
A2-80071-DF120 AX-80071-DF120	X				1.2		X	X	X	X	105"	50

- Administration set fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP
- Single-use sets, universal spike, spike protector, drip chamber, male luer-lock and tip protector
- Needleless injection sites
- Drip chamber = 20 drops/ mL

ADMINISTRATION SETS

All of the administration listed below are single-use sets for use with Zyno Medical Z-800 Infusion pumps. The administration sets (A2/AX primary administration sets, secondary administration sets, and extension administration sets) fluid paths are sterile, non-Pyrogenic, are not made with natural rubber latex, and are not made with DEHP.

A2/AX PRIMARY ADMINISTRATION SETS

PRODUCT ID	DESCRIPTION	LENGTH	# OF Y-SITES	QTY
A2-80070	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5) -Male Luer lock adaptor, (6)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	0	50
or				
AX-80070				
A2-80071-D	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-Slide clamp, (6)-Needleless injection site, (7)-Male Luer lock adaptor, (8)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	1 (distal)	50
or				
AX-80071-D				
A2-80071-P	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5)-Needleless injection site, (6)-Roller clamp, (7)-Male Luer lock adaptor, (8)-Luer lock tip protector. Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	1 (proximal)	50
or				
AX-80071-P				

A2-80072



or

(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5 & 8)-Needleless injection sites, (6)-Roller clamp, (7)-Slide clamp, (9)-Male Luer lock adaptor, (10)-Luer lock tip protector.

Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL

105"
(266cm)






2

50



AX-80072





A2/AX PRIMARY ADMINISTRATION SETS WITH 0.2 MICRON FILTER

PRODUCT ID	DESCRIPTION	LENGTH	# OF Y-SITES	QTY
A2-70071-DF	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-0.2 micron filter, (6)-Slide clamp, (7)-Needleless injection site, (8)-Male Luer lock adaptor, (9)-Luer lock tip protector.</p> <p>Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	1 (distal)	50
or				
AX-70071-DF				
A2-80071-PF	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5)-Needleless injection site, (6)-Roller clamp, (7)-0.2 micron filter, (8)-Male Luer lock adaptor, (9)-Luer lock tip protector.</p> <p>Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	1 (proximal)	50
A2-80072-F	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Back check valve, (5 & 9)-Needleless injection sites, (6)-Roller clamp, (7)-0.2 micron filter, (8)-Slide clamp, (10)-Male Luer lock adaptor, (11)-Luer lock tip protector.</p> <p>Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>	105" (266cm)	2	50
or				
AX-80072-F				

A2/AX PRIMARY ADMINISTRATION SETS WITH 1.2 MICRON FILTER

PRODUCT ID	DESCRIPTION	LENGTH	# OF Y-SITES	QTY
A2-80071-DF120	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-1.2 micron filter, (6)-Slide clamp, (7)-Needleless injection site, (8)-Male Luer lock adaptor, (9)-Luer lock tip protector.</p>	105" (266cm)	1 (distal)	50
or	<p>Priming Volume: 18 mL (approx.)/ Drop Size: 20 drops/mL</p>			
AX-80071-DF120				

SECONDARY ADMINISTRATION SETS

PRODUCT ID	DESCRIPTION	LENGTH	# OF Y-SITES	QTY
AA-80075	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-Male Luer lock adaptor, (6)-Luer lock tip protector, (7)-Bag hanger (11 inch).</p>	40" (101cm)	0	50
AX-80075	 <p>(1)-Spike tip protector, (2)-Universal spike, (3)-Drip chamber, (4)-Roller clamp, (5)-Male Luer lock adaptor, (6)-Luer lock tip protector, (7)-Bag hanger (16.5 inch).</p>	40" (101cm)	0	50



**Ready to Order? Need Help?
Let's get in touch.**

OUR CUSTOMER SERVICE TEAM

P 1 (508) 650 2008

F 1 (508) 650-2006

E orders@zynosolutions.com

www.zynosolutions.com

Hours: 8 AM - 5 PM, Monday to Friday

For technical support, please contact 866-395-1988 or techsupport@zynosolutions.com.



P +1 508 650 2008

E info@zynosolutions.com

A 177 Pine Street Natick MA 01760 USA
www.zynosolutions.com

ZYS-CTAX-002 Rev B